

Creatinine (fluid) (FCR)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

27 July 2026 20:11



TEST OVERVIEW

Test Name	Creatinine (fluid)
Test Code	FCR
Short Description	Creatinine (fluid)

OVERVIEW

Test Name	Creatinine (fluid)
Test Code	FCR
Category	Fluid Biochemistry
TAT	Main Lab: 6 Hour(s) Family Site: <6hrs
Specimen(s)	1 x Fluid - 5 mL Steril container - Red - Fluid

SPECIMEN(S)

Fluid

Specimen Type	Fluid
Specimen Format	Steril container
Specimen Colour	Red
Specimen Volume	5 mL
Sampling Order	0
Origin	Fluid
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	-
TAT	Main Lab: 6 Hour(s)

Test Stability

Family Site: <6hrs

Room Temp: 7 Day(s)
2–8°C: 7 Day(s)**CLINICAL INFORMATION****Creatinine (fluid)**

Methodology	-
Specimen Type	Fluid
Delay before pre-treatment	-
Transport Temperature	15-25°C
Transport Stability at room temp	7 Day
Transport Stability at 2–8°C	7 Day
Haemolysis interference	No

PATIENT INFORMATION

Clinical Information Required	-
Patient Collection Notes	-

COMMENTS & NOTES

LOINC Code	386-4, 25386-4
Outwork	Yes