

Erythrocyte Sedimentation Rate (ESR) (ESR)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

04 May 2026 22:49



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Erythrocyte Sedimentation Rate (ESR)
Test Code	ESR
Short Description	ESR

OVERVIEW

Test Name	Erythrocyte Sedimentation Rate (ESR)
Test Code	ESR
Category	Haematology
TAT	Main Lab: 6 Hour(s) Family Site: <8hrs
Specimen(s)	1 x Venous blood - 5 mL Tube - Lavender - EDTA Whole Blood 1 x - - 5 mL Tube - Black - ESR 1 x - - - mL Specific tube - Black - HUMASED

SPECIMEN(S)

EDTA Whole Blood

Specimen Type	EDTA Whole Blood
Specimen Format	Tube
Specimen Colour	Lavender
Specimen Volume	5 mL
Sampling Order	4
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	-

TAT Main Lab: 6 Hour(s)
Family Site: <8hrs

Test Stability Room Temp: 6 Hour(s)
2–8°C: 24 Hour(s)

ESR

Specimen Type ESR

Specimen Format Tube

Specimen Colour Black

Specimen Volume 5 mL

Sampling Order 7

Origin -

Collection time after baseline -

Transport Temperature 15-25°C

Accepted Other Specimens -

TAT Main Lab: 6 Hour(s)
Family Site: <8hrs

Test Stability Room Temp: 6 Hour(s)
2–8°C: 24 Hour(s)

HUMASED

Specimen Type HUMASED

Specimen Format Specific tube

Specimen Colour Black

Specimen Volume - mL

Sampling Order 7

Origin -

Collection time after baseline -

Transport Temperature 15-25°C

Accepted Other Specimens -

TAT Main Lab: 6 Hour(s)
Family Site: <8hrs

Test Stability Room Temp: 6 Hour(s)
2–8°C: 24 Hour(s)

PATIENT INFORMATION

Clinical Information Required -

Patient Collection Notes -

COMMENTS & NOTES

LOINC Code 37-7, 4537-7, 82477-1, 30341-2

Outwork No