

Food Screen Fx5 (Pediatrics) (FX5PAED)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

27 July 2026 20:10



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Food Screen Fx5 (Pediatrics)
Test Code	FX5PAED
Short Description	Food Screen Fx5 (Pediatrics)

OVERVIEW

Test Name	Food Screen Fx5 (Pediatrics)
Test Code	FX5PAED
Category	Allergy/Autoimmune
TAT	Main Lab: 1 Day(s) Family Site: 1 Day
Specimen(s)	1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube

SPECIMEN(S)

SST-Serum Separator Tube

Specimen Type	SST-Serum Separator Tube
Specimen Format	Tube
Specimen Colour	Gold
Specimen Volume	5 mL
Sampling Order	2
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	Serum
TAT	Main Lab: 1 Day(s)

Test Stability

Family Site: 1 Day

Room Temp: 24 Hour(s)
2–8°C: 24 Hour(s)**CLINICAL INFORMATION****Food Screen Fx5 (Pediatrics)**

Methodology	-
Specimen Type	SST-Serum Separator Tube
Delay before pre-treatment	6
Transport Temperature	15-25°C
Transport Stability at room temp	24 Hours
Transport Stability at 2–8°C	24 Hours
Haemolysis interference	No

PATIENT INFORMATION

Clinical Information Required	-
Patient Collection Notes	-

COMMENTS & NOTES

LOINC Code	-
Outwork	Yes