

Glucose 2 Hours (GLU2HR)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

04 July 2026 13:02



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Glucose 2 Hours
Test Code	GLU2HR
Short Description	Glucose 2 Hours

OVERVIEW

Test Name	Glucose 2 Hours
Test Code	GLU2HR
Category	Biochemistry
TAT	Main Lab: 2 Hour(s) Family Site: <3hrs, <2hrs
Specimen(s)	1 x Venous blood - 5 mL Tube - Grey - Sodium Fluorid Plasma 2Hrs

SPECIMEN(S)

Sodium Fluorid Plasma 2Hrs

Specimen Type	Sodium Fluorid Plasma 2Hrs
Specimen Format	Tube
Specimen Colour	Grey
Specimen Volume	5 mL
Sampling Order	5
Origin	Venous blood
Collection time after baseline	2
Transport Temperature	2-8°C
Accepted Other Specimens	-
TAT	Main Lab: 2 Hour(s)

Test Stability

Family Site: <3hrs, <2hrs

Room Temp: 2 Day(s)
2-8°C: 7 Day(s)**PATIENT INFORMATION****Clinical Information Required** -**Patient Collection Notes** -**COMMENTS & NOTES****LOINC Code** 759-5, 14759-5**Outwork**