

## Haemoglobin (HB)

Source:

Cerba Lancet Africa

View our website:

[cerbalancetafrica.com](http://cerbalancetafrica.com)

Date & Time of Export:

15 April 2026 06:08



[Click here to view the results online](#)

### TEST OVERVIEW

|                          |             |
|--------------------------|-------------|
| <b>Test Name</b>         | Haemoglobin |
| <b>Test Code</b>         | HB          |
| <b>Short Description</b> | HB          |

### OVERVIEW

|                    |  |
|--------------------|--|
| <b>Test Name</b>   | Haemoglobin  |
| <b>Test Code</b>   | HB   |
| <b>Category</b>    | Haematology  |
| <b>TAT</b>         | Main Lab: 2 Hour(s)<br>Family Site: <2hrs                  |
| <b>Specimen(s)</b> | 1 x Venous blood - 5 mL Tube - Lavender - EDTA Whole Blood |

### SPECIMEN(S)

#### EDTA Whole Blood

|                                       |                     |
|---------------------------------------|---------------------|
| <b>Specimen Type</b>                  | EDTA Whole Blood    |
| <b>Specimen Format</b>                | Tube                |
| <b>Specimen Colour</b>                | Lavender            |
| <b>Specimen Volume</b>                | 5 mL                |
| <b>Sampling Order</b>                 | 4                   |
| <b>Origin</b>                         | Venous blood        |
| <b>Collection time after baseline</b> | -                   |
| <b>Transport Temperature</b>          | 15-25°C             |
| <b>Accepted Other Specimens</b>       | -                   |
| <b>TAT</b>                            | Main Lab: 2 Hour(s) |

**Test Stability**

Family Site: &lt;2hrs

Room Temp: 24 Hour(s)  
2-8°C: 24 Hour(s)**PATIENT INFORMATION****Clinical Information Required** -**Patient Collection Notes** -**COMMENTS & NOTES****LOINC Code** 8-7, 718-7**Outwork**