

HIV Rapid Screen (HIVR)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

27 July 2026 20:16



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	HIV Rapid Screen
Test Code	HIVR
Short Description	HIV Rapid

OVERVIEW

Test Name	HIV Rapid Screen
Test Code	HIVR
Category	Immunoassay
TAT	Main Lab: 2 Hour(s) Family Site: <3hrs
Specimen(s)	1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube

SPECIMEN(S)

SST-Serum Separator Tube

Specimen Type	SST-Serum Separator Tube
Specimen Format	Tube
Specimen Colour	Gold
Specimen Volume	5 mL
Sampling Order	2
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	2-8°C
Accepted Other Specimens	EDTA Whole Blood EDTA Plasma Serum

TAT

Main Lab: 2 Hour(s)
Family Site: <3hrs

Test Stability

Room Temp: 1 Day(s)
2–8°C: 8 Day(s)

PATIENT INFORMATION

Clinical Information Required

-

Patient Collection Notes

-

COMMENTS & NOTES

LOINC Code

580-4, 49580-4

Outwork

No