

Iron (IRON)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

15 April 2026 06:07



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Iron
Test Code	IRON
Short Description	Iron

OVERVIEW

Test Name	Iron
Test Code	IRON
Category	Biochemistry
TAT	Main Lab: 6 Hour(s) Family Site: <8hrs, <6hrs
Specimen(s)	1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube

SPECIMEN(S)

SST-Serum Separator Tube

Specimen Type	SST-Serum Separator Tube
Specimen Format	Tube
Specimen Colour	Gold
Specimen Volume	5 mL
Sampling Order	2
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	Lithium Heparin Plasma Sodium Heparin Plasma Serum

TAT

Main Lab: 6 Hour(s)
Family Site: <8hrs, <6hrs

Test Stability

Room Temp: 10 Hour(s)
2–8°C: 7 Day(s)

PATIENT INFORMATION

Clinical Information Required -

Patient Collection Notes -

COMMENTS & NOTES

LOINC Code 798-3, 14798-3

Outwork