

Magnesium (MAG)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

15 April 2026 06:01



[Click here to view the results online](#)

TEST OVERVIEW

| | |
|--------------------------|-----------|
| Test Name | Magnesium |
| Test Code | MAG |
| Short Description | MG |

OVERVIEW

| | |
|--------------------|--|
| Test Name | Magnesium |
| Test Code | MAG |
| Category | Biochemistry |
| TAT | Main Lab: 6 Hour(s) Family Site: <8hrs, <6hrs |
| Specimen(s) | 1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube |

SPECIMEN(S)

SST-Serum Separator Tube

| | |
|---------------------------------------|--|
| Specimen Type | SST-Serum Separator Tube |
| Specimen Format | Tube |
| Specimen Colour | Gold |
| Specimen Volume | 5 mL |
| Sampling Order | 2 |
| Origin | Venous blood |
| Collection time after baseline | - |
| Transport Temperature | 15-25°C |
| Accepted Other Specimens | Serum Lithium Heparin Plasma Sodium Heparin Plasma |

TAT

Main Lab: 6 Hour(s)
Family Site: <8hrs, <6hrs

Test Stability

Room Temp: 7 Day(s)
2–8°C: 7 Day(s)

PATIENT INFORMATION

Clinical Information Required -

Patient Collection Notes -

COMMENTS & NOTES

LOINC Code 01-3, 2601-3

Outwork