

Malaria (Thin Smear) (MALT)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

01 April 2026 21:26



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Malaria (Thin Smear)
Test Code	MALT
Short Description	Malaria (Thin Smear)

OVERVIEW

Test Name	Malaria (Thin Smear)
Test Code	MALT
Category	Haematology
TAT	Main Lab: 12 Hour(s) Family Site: <12hrs
Specimen(s)	1 x Venous blood - 5 mL Tube - Lavender - EDTA Whole Blood

SPECIMEN(S)

EDTA Whole Blood

Specimen Type	EDTA Whole Blood
Specimen Format	Tube
Specimen Colour	Lavender
Specimen Volume	5 mL
Sampling Order	4
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	Capillary Blood
TAT	Main Lab: 12 Hour(s)

Test Stability

Family Site: <12hrs
Room Temp: 12 Hour(s)
2-8°C: 24 Hour(s)

PATIENT INFORMATION**Clinical Information Required**

Visited country
Infection signs
Start date of clinical signs
Use of malaria prophylaxis AND/OR treatment

Patient Collection Notes

-

COMMENTS & NOTES**LOINC Code**

271-8, 33271-8

Outwork

No