

Cardiac Troponin (Troponin)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

19 January 2026 06:30



TEST OVERVIEW

Test Name	Cardiac Troponin
Test Code	Troponin
Short Description	Troponin

OVERVIEW

Test Name	Cardiac Troponin
Test Code	Troponin
Category	POC Cardiac
TAT	Main Lab: 6 Hour(s) Family Site: <8hrs
Specimen(s)	1 x Venous blood - 5 mL Tube - Green - Lithium Heparin Whole Blood 1 x - - 5 mL Tube - Gold - SST-Serum Separator Tube 1 x - - 5 mL Tube - Lavender - EDTA Plasma

SPECIMEN(S)

Lithium Heparin Whole Blood

Specimen Type	Lithium Heparin Whole Blood
Specimen Format	Tube
Specimen Colour	Green
Specimen Volume	5 mL
Sampling Order	3
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	Serum

	Lithium Heparin Plasma EDTA Plasma Sodium Heparin Plasma SST-Serum Separator Tube
TAT	Main Lab: 6 Hour(s) Family Site: <8hrs

Test Stability	Room Temp: 8 Hour(s) 2–8°C: -
-----------------------	----------------------------------

SST-Serum Separator Tube

Specimen Type	SST-Serum Separator Tube
Specimen Format	Tube
Specimen Colour	Gold
Specimen Volume	5 mL
Sampling Order	2
Origin	-
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	Serum Lithium Heparin Plasma EDTA Plasma Sodium Heparin Plasma SST-Serum Separator Tube

TAT	Main Lab: 6 Hour(s) Family Site: <8hrs
------------	---

Test Stability	Room Temp: 8 Hour(s) 2–8°C: -
-----------------------	----------------------------------

EDTA Plasma

Specimen Type	EDTA Plasma
Specimen Format	Tube
Specimen Colour	Lavender
Specimen Volume	5 mL
Sampling Order	4
Origin	-
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	Serum Lithium Heparin Plasma EDTA Plasma Sodium Heparin Plasma SST-Serum Separator Tube

TAT	Main Lab: 6 Hour(s) Family Site: <8hrs
------------	---

Test Stability	Room Temp: 8 Hour(s) 2–8°C: -
-----------------------	----------------------------------

PATIENT INFORMATION

Clinical Information Required -

Patient Collection Notes -

COMMENTS & NOTES

LOINC Code 97-9, 6597-9, 10839-9, 89579-7

Outwork