

Toxoplasma IgM (TOXM)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

15 April 2026 06:08



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Toxoplasma IgM
Test Code	TOXM
Short Description	TOX IgM

OVERVIEW

Test Name	Toxoplasma IgM
Test Code	TOXM
Category	Serology
TAT	Main Lab: 24 Hour(s) Family Site: 1 Day(s)
Specimen(s)	1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube

SPECIMEN(S)

SST-Serum Separator Tube

Specimen Type	SST-Serum Separator Tube
Specimen Format	Tube
Specimen Colour	Gold
Specimen Volume	5 mL
Sampling Order	2
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	2-8°C
Accepted Other Specimens	Lithium Heparin Plasma Serum EDTA Plasma

	Sodium Heparin Plasma
TAT	Main Lab: 24 Hour(s) Family Site: 1 Day(s)
Test Stability	Room Temp: 3 Day(s) 2–8°C: 14 Day(s)

PATIENT INFORMATION

Clinical Information Required	Date of pregnancy Pregnancy
Patient Collection Notes	-

COMMENTS & NOTES

LOINC Code	262-2, 12262-2, 25542-2
Outwork	<input type="button" value="No"/>