

Total Protein (TP)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

04 July 2026 12:48



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Total Protein
Test Code	TP
Short Description	TP

OVERVIEW

Test Name	Total Protein
Test Code	TP
Category	Biochemistry
TAT	Main Lab: 4 Hour(s) Family Site: <5hrs, <4hrs
Specimen(s)	1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube

SPECIMEN(S)

SST-Serum Separator Tube

Specimen Type	SST-Serum Separator Tube
Specimen Format	Tube
Specimen Colour	Gold
Specimen Volume	5 mL
Sampling Order	2
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	Lithium Heparin Plasma Serum Sodium Heparin Plasma

TAT

EDTA Plasma

Main Lab: 4 Hour(s)
Family Site: <5hrs, <4hrs

Test Stability

Room Temp: 1 Day(s)
2–8°C: 1 Month(s)

PATIENT INFORMATION

Clinical Information Required -

Patient Collection Notes -

COMMENTS & NOTES

LOINC Code 85-2, 2885-2

Outwork