

Adrenocortotropin Hormone (ACTH)-Midnight (ACTH-MIDNIGHT)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

27 July 2026 20:15



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Adrenocortotropin Hormone (ACTH)-Midnight
Test Code	ACTH-MIDNIGHT
Short Description	Adrenocortotropin Hormone (ACTH)-Midnight

OVERVIEW

Test Name	Adrenocortotropin Hormone (ACTH)-Midnight
Test Code	ACTH-MIDNIGHT
Category	Immunoassay
TAT	Main Lab: 12 Hour(s) Family Site: <12hrs
Specimen(s)	1 x Venous blood - 5 mL Tube - Lavender - EDTA Plasma

SPECIMEN(S)

EDTA Plasma

Specimen Type	EDTA Plasma
Specimen Format	Tube
Specimen Colour	Lavender
Specimen Volume	5 mL
Sampling Order	4
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	-80°C
Accepted Other Specimens	-

TAT

Main Lab: 12 Hour(s)
Family Site: <12hrs

Test Stability

Room Temp: 2 Hour(s)
2–8°C: 5 Hour(s)

PATIENT INFORMATION

Clinical Information Required

-

Patient Collection Notes

Collection around midnight

Immediately after collection, place the tube in an ice bath.

COMMENTS & NOTES

LOINC Code

092-1, 48092-1

Outwork

No