

## Calcium (24 Hrs Urine) (CA24)

Source:

Cerba Lancet Africa

View our website:

[cerbalancetafrica.com](http://cerbalancetafrica.com)

Date & Time of Export:

15 April 2026 06:00



[Click here to view the results online](#)

### TEST OVERVIEW

<b>Test Name</b>	Calcium (24 Hrs Urine)
<b>Test Code</b>	CA24
<b>Short Description</b>	Calcium (24 Hrs Urine)

### OVERVIEW

<b>Test Name</b>	Calcium (24 Hrs Urine)
<b>Test Code</b>	CA24
<b>Category</b>	Urine Biochemistry
<b>TAT</b>	Main Lab: 4 Hour(s) Family Site: <4hrs
<b>Specimen(s)</b>	1 x Urine - 2000 mL 24-Hours urine container - Red - Urine 24hrs - HCL

### SPECIMEN(S)

#### Urine 24hrs - HCL

<b>Specimen Type</b>	Urine 24hrs - HCL
<b>Specimen Format</b>	24-Hours urine container
<b>Specimen Colour</b>	Red
<b>Specimen Volume</b>	2000 mL
<b>Sampling Order</b>	0
<b>Origin</b>	Urine
<b>Collection time after baseline</b>	-
<b>Transport Temperature</b>	-
<b>Accepted Other Specimens</b>	-
<b>TAT</b>	Main Lab: 4 Hour(s)

**Test Stability**

Family Site: &lt;4hrs

Room Temp: -  
2-8°C: -**PATIENT INFORMATION****Clinical Information Required**24-hour urine volume (mL)  
Collection start date  
Collection start time  
Collection end time  
Collection End date**Patient Collection Notes**

24-hour urine collection:

- On the first day, on waking, eliminate the first urine in the toilet (note the date and time).
- For the next 24 hours, collect all the urine from the day and night, including that from the next morning when you get up at the same time.
- Recap and store the container between 2 and 8 degrees C. between each micturition.
- Note the date and time of the end of micturition on the bottle.
- Be sure to close the container securely for transport.
- Bring all the urine to the laboratory as soon as possible after collection.

**COMMENTS & NOTES****LOINC Code**

362-5, 25362-5

**Outwork****No**