

aPTT-Correction (APTT CORR)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

27 July 2026 20:11



TEST OVERVIEW

Test Name	aPTT-Correction
Test Code	APTT CORR
Short Description	aPTT-Correction

OVERVIEW

Test Name	aPTT-Correction
Test Code	APTT CORR
Category	Haemostasis
TAT	Main Lab: 3, 4 Hour(s) Family Site: <5hrs, <4hrs
Specimen(s)	1 x Venous blood - 5 mL Tube - Light Blue - Citrated Plasma

SPECIMEN(S)

Citrated Plasma

Specimen Type	Citrated Plasma
Specimen Format	Tube
Specimen Colour	Light Blue
Specimen Volume	5 mL
Sampling Order	1
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	-
Accepted Other Specimens	-
TAT	Main Lab: 3, 4 Hour(s)

Test Stability

Family Site: <5hrs, <4hrs

Room Temp: -
2–8°C: -**CLINICAL INFORMATION****Rosner Index**

Methodology	Calculation
Specimen Type	Citrated Plasma
Delay before pre-treatment	-
Transport Temperature	-
Transport Stability at room temp	- -
Transport Stability at 2–8°C	- -
Haemolysis interference	No

APTT (Patient)

Methodology	Clotting: Waterbath method
Specimen Type	Citrated Plasma
Delay before pre-treatment	-
Transport Temperature	-
Transport Stability at room temp	- -
Transport Stability at 2–8°C	- -
Haemolysis interference	No

APTT (Control)

Methodology	Clotting: Waterbath method
Specimen Type	Citrated Plasma
Delay before pre-treatment	-
Transport Temperature	-
Transport Stability at room temp	- -
Transport Stability at 2–8°C	- -
Haemolysis interference	No

aPTT-Mix (patient+normal plasma)

Methodology	Clotting: Waterbath method
Specimen Type	Citrated Plasma

Delay before pre-treatment -

Transport Temperature -

Transport Stability at room temp - -

Transport Stability at 2–8°C - -

Haemolysis interference

PATIENT INFORMATION

Clinical Information Required -

Patient Collection Notes -

COMMENTS & NOTES

LOINC Code 46-9, 5946-9,

Outwork