

Endomysium Anti IgA (ENDOMY-ANTI-IgA)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

15 April 2026 06:08



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Endomysium Anti IgA
Test Code	ENDOMY-ANTI-IgA
Short Description	Endomy-Anti-IgA

OVERVIEW

Test Name	Endomysium Anti IgA
Test Code	ENDOMY-ANTI-IgA
Category	Referred tests
TAT	Main Lab: 48 Hour(s) Family Site: -
Specimen(s)	1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube

SPECIMEN(S)

SST-Serum Separator Tube

Specimen Type	SST-Serum Separator Tube
Specimen Format	Tube
Specimen Colour	Gold
Specimen Volume	5 mL
Sampling Order	2
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	-
Accepted Other Specimens	Serum
TAT	Main Lab: 48 Hour(s)

Test Stability

Family Site: -

Room Temp: -
2-8°C: -

PATIENT INFORMATION

Clinical Information Required -

Patient Collection Notes -

COMMENTS & NOTES

LOINC Code 863-9, 10863-9

Outwork Yes