

## Aspergillus Galactomannan Antigen (ASPERG GALAC)

Source:

Cerba Lancet Africa

View our website:

[cerbalancetafrica.com](http://cerbalancetafrica.com)

Date & Time of Export:

15 April 2026 06:03



[Click here to view the results online](#)

### TEST OVERVIEW

<b>Test Name</b>	Aspergillus Galactomannan Antigen
<b>Test Code</b>	ASPERG GALAC
<b>Short Description</b>	Asperg Galac

### OVERVIEW

<b>Test Name</b>	Aspergillus Galactomannan Antigen
<b>Test Code</b>	ASPERG GALAC
<b>Category</b>	Referred tests
<b>TAT</b>	Main Lab: 7 Day(s) Family Site: -
<b>Specimen(s)</b>	1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube

### SPECIMEN(S)

#### SST-Serum Separator Tube

<b>Specimen Type</b>	SST-Serum Separator Tube
<b>Specimen Format</b>	Tube
<b>Specimen Colour</b>	Gold
<b>Specimen Volume</b>	5 mL
<b>Sampling Order</b>	2
<b>Origin</b>	Venous blood
<b>Collection time after baseline</b>	-
<b>Transport Temperature</b>	-
<b>Accepted Other Specimens</b>	Serum
<b>TAT</b>	Main Lab: 7 Day(s)

**Test Stability**

Family Site: -

Room Temp: -  
2-8°C: -

### PATIENT INFORMATION

**Clinical Information Required** -

**Patient Collection Notes** -

### COMMENTS & NOTES

**LOINC Code** 383-9, 35383-9

**Outwork**  Yes