

Gastric Antibodies (Anti-Parietal) (GASTRIC_AB)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

27 July 2026 20:13



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Gastric Antibodies (Anti-Parietal)
Test Code	GASTRIC_AB
Short Description	Gastric Antibodies (Anti-Parietal)

OVERVIEW

Test Name	Gastric Antibodies (Anti-Parietal)
Test Code	GASTRIC_AB
Category	Allergy/Autoimmune
TAT	Main Lab: 72 Hour(s) Family Site: 3 Days
Specimen(s)	1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube

SPECIMEN(S)

SST-Serum Separator Tube

Specimen Type	SST-Serum Separator Tube
Specimen Format	Tube
Specimen Colour	Gold
Specimen Volume	5 mL
Sampling Order	2
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	Serum
TAT	Main Lab: 72 Hour(s)

Test Stability

Family Site: 3 Days

Room Temp: 4 Day(s)
2–8°C: 7 Day(s)**PATIENT INFORMATION****Clinical Information Required**

-

Patient Collection Notes

-

COMMENTS & NOTES**LOINC Code**

969-6, 26969-6

Outwork

Yes