

Estradiol (E2)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

15 April 2026 06:07



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Estradiol
Test Code	E2
Short Description	E2

OVERVIEW

Test Name	Estradiol
Test Code	E2
Category	Immunoassay
TAT	Main Lab: 4 Hour(s) Family Site: <4hrs
Specimen(s)	1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube

SPECIMEN(S)

SST-Serum Separator Tube

Specimen Type	SST-Serum Separator Tube
Specimen Format	Tube
Specimen Colour	Gold
Specimen Volume	5 mL
Sampling Order	2
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	Serum Lithium Heparin Plasma EDTA Plasma

TAT	Sodium Heparin Plasma Main Lab: 4 Hour(s) Family Site: <4hrs
Test Stability	Room Temp: 8 Hour(s) 2–8°C: 24 Hour(s)

PATIENT INFORMATION

Clinical Information Required	- Day of menstrual cycle il vitro fertilization
Patient Collection Notes	-

COMMENTS & NOTES

LOINC Code	097-6, 83097-6
Outwork	<input type="button" value="No"/>